

Finishing Accessories For Regal Cage

Thermometer/hygrometer monitors temperature and humidity.



Melted ice used for cooling oxygen is drained from compartment through convenient valve.

Clear polycarbonate panel resists scratches and discoloration.

Therapy Doors

Transforms Regal Cage into a respiratory therapy cage. Allows for the following therapies:

- 1) Oxygen
- 2) Oxygen with moisture
- 3) Oxygen with moisture and medication
- 4) Compressed air with moisture
- 5) Compressed air with moisture and medication

Clear door allows for patient observation and features rubber gasket seal which is open at the corners to create a positive atmosphere, yet allows for venting. A continuous-flow nebulizer mounted on door delivers a regulated flow of the preferred therapy. Ice compartment mounted on door holds ice to cool therapy liquid, if needed. Thermometer/ hygrometer, mounted on door interior, indicates temperature and humidity of cage. Door mounts hinged left onto front of cage in place of the regular door; uses regular hinges and latches to close securely with .25-in. (.64cm) inner diameter x .062-in. (.16cm) thick tubing. Only available for single-door cages.

| P/N | FITS CAGE SIZE | |
|-------------------------|-----------------------------|------------------|
| | Inches | Centimeters |
| 12155-00-DRDRAA | 24W x 24H | 60.96W x 60.96H |
| 12155-00-DREPAA | 24W x 30H | 60.96W x 76.20H |
| 12155-00-EPEPAA | 30W x 30H | 76.20W x 76.20H |
| 12155-00-FNEPAA | 36W x 30H | 91.44W x 76.20H |
| 12155-00-HJEPAA | 48W x 30H | 121.92W x 76.20H |
| Additional Parts | | |
| 212126 | Nebulizers with Supply Hose | |

