

Steam Distiller

Purification System, Model DS1000

Pure Clean Water, Made Simple...

Simple one touch operation.

You will love using our model DS1000 Steam Distiller. Fill the stainless-steel boiler with tap water and press the button. Shortly after you will have the purest possible water.

Features:

- Produces 1 gallon of 100% steam distilled water in 4 hours (Up to 6 gallons per day)
- Reduces 99% of total dissolved solids (TDS)
- Dripless stop-and-serve feature



PN# DS1000


Tuttnauer
Innovation · Legacy · Partnership

How It Works:

The Tuttnauer Steam Distiller Model DS1000 system produces purified water by steam distillation. As the water temperature rises to 212°F (100°C), it kills bacteria, cysts, and viruses that may be present.

1. As the boiled water
2. Converts to steam
3. It rises and leaves behind dead microbes, and other contaminants and low boiling light gases are discharged through the gaseous vent
4. Steam is then cooled in the stainless-steel condenser turning it into high purity water.
5. Distilled water percolates through a final polishing coconut shell carbon filter
6. It is then collected in the Tritan™ carafe which is BPA free.



More Features:

- Automatic shutoff after each 1 gallon cycle
- 800 watts of power, 120 volt 60 hZ, fits standard outlet
- Effectively removes arsenic, nitrates, sodium, chlorine, bacteria and viruses.
- 1-gallon (4 liter) BPA Free carafe included
- Item Weight 11.43 pounds
- Dimensions: 15.25" H x 15" W x 9.75" D



The seamless, stainless steel boiler has a concealed heating element that won't become coated with scale and mineral deposits. A simple twist locks the 7 inch stainless steel cover firmly into place.

Accessory Items:



6-Pack Certified Carbon Filter Bags
PN# 8606



Extra Tritan™ BPA-FREE Carafe
PN# 3250



Kleenwise™ Distiller Cleaner and De-Scaler
PN# KSD

Tuttnauer
Innovation · Legacy · Partnership

Learn from our experts.

Visit our blog: tuttnauer.com/blog



Tuttnauer USA Co. Ltd.

25 Power Drive, Hauppauge, NY 11788

Tel: +800 624 5836, +631 737 4850, Fax: +631 737 0720

E-mail: info@tuttnauerUSA.com, www.tuttnauerUSA.com