

# Avante NXT5

## Infusion Pump

### More precise, safer and easier to use

The Avante NXT5 contains a sophisticated motor and dual detectors to guarantee infusion precision and stability. Intelligent pulse compensation technology and heating technology ensured the consistency of long time infusion rate.

#### FEATURES

- High precision pressure sensor is able to detect down to 50mmHg within 1ms, alarm will be triggered once blocking happens.
- Dual CPU monitors the whole infusion process in real time.
- Built-in drip sensor to control fluids flux.
- 800 infusion records and last configuration memory.
- Built-in Lithium battery, over 8 hours working time after fully charged.
- Rate Mode, Drop Mode and Time Mode for wide range of applications.
- 4" LCD display with big font
- Ergonomics operation interface make it easy to adjust and monitor the parameters so as to reduce mistake operation.
- Micro infusion speed 0.1ml/h up to high speed 1500ml/h satisfies special clinic scenario including requirements for Neonatal.
- Software adjusting function suit for any brand of infusion sets.
- Automatic light up while door was opened, easy to operate at night.



#### SPECIFICATIONS



**Weight:**  
3.97 lb (1.8 kg)



**Display:**  
6.85 x 4.96 x 8.46 in  
(174 x 126 x 215 mm)



**Power Supply:**  
11.1 V DC

**Battery:**  
Built-in Lithium,  
over 8 hours on a full charge



**Infusion Modes:**  
Rate Mode, Drop Mode  
and Time Mode

**Micro Infusion Speed:**  
0.1ml/h up to 1500ml/h

**Accuracy:** < ±3%




2375-2020-12-04



Intelligent pulse compensation technology is to compensate infusion fluctuation and accurately adjust the infusion speed to realize safety and stable infusion.

The infusion tube will harden and slow to rebound when under low temperature, which lower down infusion precision in clinic application. Intelligent heating technology is able to provide adequate temperature for infusion system from 25°C~40°C with over temperature protection, alleviate discomfort of patient.

Horizontal door clip for easier close and open the door.



LED alarm lamp, infusion state can be seen from a long distance. Variety of audible and visible alarms to guarantee safety infusion.

Ultrasound probe can detect 25ul bubble to prevent it from entering into an animal and leading to aeroembolism; separate design make the tube installation easier and also reduce the risk of tube damage.

Patented anti-bolus system to reduce significantly bolus after occlusion sudden release. Automatic detect door closing situation can effectively prevent drug flowing and protect the tube.